


MCR-SEAT

Mobility Care Range Shower Seat

meir®



EXPLODED VIEW



Item No.	Part Name	Qty.
1	Seat	1
2	Bracket Back Cover	1
3	Cover Strip	1
4	Slow Down Damper	2
5	Rotating Shaft	2
6	Mounting Base	1
7	Base Cover	1
8	Nylon Bolt Cap	6
9	Gasket	6
10	Screw	6

CARE

- Do not apply physical items (such as tools) directly to the product.
- Do not install the products using any form of acetone silicones.
- Never use house-hold or commercial detergents, citrus based cleaners, or abrasive cleaners.
- Clean using a mild washing soap solution rinsed, and dry with a soft cloth.
- Refer to the installation manual for correct installation.

FEATURES

- Refer to the Meir website for warranty terms & conditions, and available finishes.
- Load capacity of 150 kg when installed correctly.
- Soft opening seat movement.

SPECIFICATIONS

