

MCR-SEAT
Mobility Care Range Shower Seat

meir®



EXPLODED VIEW



Item No.	Part Name	Qty.
1	Seat	1
2	Bracket Back Cover	1
3	Cover Strip	1
4	Slow Down Damper	2
5	Rotating Shaft	2
6	Mounting Base	1
7	Base Cover	1
8	Nylon Bolt Cap	6
9	Gasket	6
10	Screw	6


CARE


- Do not apply physical items (such as tools) directly to the product.
- Do not install the products using any form of acetone silicones.
- Never use house-hold or commercial detergents, citrus based cleaners, or abrasive cleaners.
- Clean using a mild washing soap solution rinsed, and dry with a soft cloth.
- Refer to the installation manual for correct installation.


FEATURES

- Load capacity of 150kg when installed correctly.
- Soft opening seat movement.

SPECIFICATIONS







– brought by –
elenni
Kitchen & Bathroom Collection

www.elenni.com.au
customerservice@elenni.com.au
Contact us 07 32099644
1 Bryants Rd, Loganholme QLD 4129